

# SARS-CoV-2 & Influenza A/B & RSV Antigen Kit (Colloidal Gold)




## Features

- Reliable test performance
- Rapid result answers in only 15 minutes
- Easily visible and readable result
- Simple operation procedure, no equipment required
- Suitable for large-scale rapid screening
- 4 in One combo test. With only one swab delivers 4 results
- Quick differentiation of SARS-CoV-2, Influenza A/B and RSV infection.

COVID-19, Influenza, RSV are all highly contagious respiratory infections caused by different viruses: The flu by influenza virus, COVID-19 by SARS-CoV-2 virus, and RSV by respiratory syncytial virus. They are all considered as highly contagious respiratory disease.

## Material provided

Sampling Swab



Test cartridge




Extraction Tube with buffer

# Specifications


Cat.No.	Pack Size
CG127001	1 Test/Kit
CG127025	25 Tests/Kit

## Operation



1

Collect sample



2

Mix liquid



3

Add sample then read results

## Result interpretation

Results	SARS-CoV-2 results (Left window)	Influenza A/B results (Middle window)	RSV results (Right window)
Positive	<p>SARS-CoV-2 Positive</p>	<p>Influenza A Positive</p>	<p>RSV Positive</p>
		<p>Influenza B Positive</p>	
		<p>Influenza A &amp; B Positive</p>	
Negative	<p>SARS-CoV-2 Negative</p>	<p>Influenza Negative</p>	<p>RSV Negative</p>
Invalid	<p>Invalid SARS-CoV-2 test</p>	<p>Invalid Influenza test</p>	<p>Invalid RSV test</p>

### CONTACT

**Goldsite Diagnostics Inc.**  
 No. 103C, 503C & 504D, Technology Building & # 3A & 4A, Technology Building Annex, Zhaoshang Sub-District, Nanshan District, Shenzhen, China, 518067

0086 755 2689 0807  
 sales@goldsite.com.cn  
 en.goldsite.com.cn



@ Goldsite Diagnostics